




GeneMarker HID




Quick Start Guide




GeneMarker HID Workflow


STR Human Identity Software

GeneMarkerHID supports numerous applications for Human Identity, including: Database Searches, Mixture Analysis, Kinship Analysis, and Paternity Testing.





Process Raw Data



Click the  icon or navigate to Project → Run to process your data:

Run Wizard


Template Selection


Set the template of the project

☒ Select an existing template or create one

- ☒ GlobalFile
- ☒ **Identifiler**
- ☒ Minifiler
- ☒ PowerPlex_16
- ☒ PowerPlex_18D
- ☒ PowerPlex_21
- ☒ PowerPlex_ESI_17
- ☒ PowerPlex_Fusion



Template Name:

Panel: 

Size Standard: 

Standard Color:

Analysis Type:

 Save  Delete

(1) Select the pre-saved template corresponding to your kit or chemistry.

(2) Or, manually select your Panel, Size Standard, and Standard color.

(3) After making your selections, click Next >>


Use default settings or customize the Data Processing Options.

Minimum Intensity and Percentage Global Max are applied to peaks called *outside* any marker ranges. Marker specific parameters may be adjusted using the Panel Editor tool (discussed later).

Run Wizard

Data Process - HID Analysis
Set data process options

Raw Data Analysis

☒ Auto Range (frame) 

Start: End:

☒ Smooth ☐ Enhanced Smooth

☐ Peak Saturation ☒ Baseline Subtraction

☐ Enhanced Baseline Subtraction

☒ Pull-up Correction ☒ Spike Removal


Size Call

☒ Local Southern ☐ Cubic Spline

Allele Call


☒ Auto Range (bps)

Start: End:

Peak Detection Threshold: 

Min Intensity: Max Intensity:

Percentage > Global Max



Run Wizard

Additional Settings - HID Analysis
Set additional options related to the different analysis type

Allelic Ladder:

Positive Control Template:

☐ Allele Evaluation

Peak Score:
Reject < Check < Pass

☐ Mixture Evaluation

Valid Mixture Peak Percentage: %

Min Mixture Marker Number:


☒ Auto Select Best Ladder


☒ Auto Panel Adjustment

<< Back Cancel

Select these two settings to automatically calibrate the ladder to the instrument run, and to enable the program to select the best fit ladder for each sample.

After making your selections click OK





Panel Editor

File Tools Help

Marker Specific settings can be set in the panel editor (Tools → Panel Editor).

Select your panel from the file tree to the left.

Right click on a marker and select **Edit Marker to view analysis thresholds for that locus.**

Select the panel from the **Panel Template list to make and save any changes to the marker parameters.**

Edit Marker

Marker Parameters:

Marker Name: D8S1179

Nucleotide Repeats (x): 4

Boundary: 118.3 To 174.0

Min Homozygote Intensity: 50

< = Inconclusive < = 200

Min Heterozygote Intensity: 50

< = Inconclusive < = 100

Max Heterozygote Imbalance: 60 %

Min Heterozygote Imbalance: 20 %

☐ Apply Homo/Hetero Settings to All Markers

Stutter Filter: N - x 8.2 %

N - 2x 1 %


N + x 1 %


☐ Apply Stutter Settings to All Markers


OK Cancel

[Panel Name]: Identifiler

[Ploidy]: 2







[Help](#)


Report Table Actions:

- Insert an Allele
- Delete an Allele
- Delete Alleles
- Edit an Allele
- Confirm an Allele
- Confirm Alleles
- Sort Alleles

Peak Quality Reasons:

- LS = Low Score
- OL = Off Ladder
- OB = Out of Bin
- BC = Bin Conflict
- SR = Saturated (Repaired)
- SP = Saturated (Pull-up)
- PL = Beyond Ploidy
- LO = Low Intensity
- HI = High Intensity
- IMB = Heterozygote Imbalance
- IHO = In Homozygote Inconclusive
- IHE = In Heterozygote Inconclusive

Quality Reason definitions can be displayed by clicking the *Help* button.




Click the  icon to generate print reports.

Enter a template name and click **Save**, below, to save settings for future reports.

Use the tabs to toggle between Standard, Advanced, and Page settings.

Click **Preview** to review the image prior to printing.



The image shows a 'Print Report' dialog box with several sections and tabs. The 'Standard' tab is selected. The 'Template Name' field contains 'Grouped by Dye'. The 'Print Type' section has 'Normal' selected. The 'Samples' section has 'All Samples' selected. The 'Contents' section has 'Electropherogram', 'PeakTable', and 'Follow Trace Chart' selected. The 'Dyes' section has 'FAM', 'VIC', 'NED', and 'PET' selected. The 'Mix Dyes' and 'Hide Bins' options are unchecked. The 'Preview' button is highlighted with a yellow arrow. The background shows a blue electropherogram trace with a peak labeled '0.97'.

Print Report

☒ Select an Existing Template ☐ Create One

Template Name:

☒ Grouped by Dye
☒ LIMS
☒ Standard

☐ Use Last Template





Standard | Advanced | Page

Print Type
☒ Normal ☐ Chart Overlay

Samples
☒ All Samples
☐ Selected Samples

Contents
☒ Electropherogram
☒ PeakTable
☒ Follow Trace Chart
☐ Start after All Charts Finished
☐ Start on Separate Page
☐ Forensics Table

Dyes
☒ FAM ☒ LIZ
☒ VIC ☐ Purple
☒ NED
☒ PET
☐ Mix Dyes
☐ Hide Bins

 Preview  Save  Delete  Ok

Reports can be printed directly, or saved as a PDF, JPEG, or PNG image.

Export Report to Files

Export Format:

PDF file

File Naming Method

☐ Named by sample name

☐ Start by Page Number

☐ Named by page number

☐ Save Group Samples as One File

Export Directory:

C:\Program Files (x86)\SoftGenetics\GeneMarker_I ...

Ok

Cancel

SoftGenetics

GeneMarker HID V2.7.1

Allele Report

2/9/2015 9:29:37 AM

Page 1

Allele Report: City_Crime_Lab

Project: Ladders OK.SGF

User: Admin


Template: Identifier

Panel: Identifier


Project Comments:

Sample 1: Ref Ladder: PAT_Ladder_3.fsa Run date and time: 02/27/2007 - 08:05:53 -> 02/27/2007 - 08:41:31


Dye: FAM - 8 peaks - PAT_10_C3.fsa




Dye: VIC - 10 peaks - PAT_10_C3.fsa





Dye: NED - 7 peaks - PAT_10_C3.fsa





































































Dye: PET - 4 peaks - PAT_10_C3.fsa



The Report Table can be saved as an excel or text file by simply clicking the  icon.

The Report Table will be saved in the same format in which it is being displayed. Click the  icon to modify the contents and format of the report table.

Report  		Bin				 Help	
	Sample	Marker	Allele#	Allele#	Allele#		
1	PAT_10_C1.fsa	D8S1179	 11	 13			
2	PAT_10_C1.fsa	D21S11	 28	 32.2			
3	PAT_10_C1.fsa	D7S820	 10	 12			
4	PAT_10_C1.fsa	CSF1PO	 11	 12			
	PAT_10_C2.fsa	D8S1179	 11	 13			
	PAT_10_C2.fsa	D21S11	 27	 32.2			
	PAT_10_C2.fsa	D7S820	 9	 10			
	PAT_10_C2.fsa	CSF1PO	 11	 13			
	PAT_10_C3.fsa	D8S1179	 12	 13			
	PAT_10_C3.fsa	D21S11	 27	 32.2			
11	PAT_10_C3.fsa	D7S820	 8	 9			
2	PAT_10_C3.fsa	CSF1PO	 11	 12			
3	PAT_10_F.fsa	D8S1179	 11	 12			
4	PAT_10_F.fsa	D21S11	 27	 28			
5	PAT_10_F.fsa	D7S820	 9	 12			
6	PAT_10_F.fsa	CSF1PO	 11	 11			
7	PAT_11_C.fsa	D8S1179	 13	 16			
8	PAT_11_C.fsa	D21S11	 24.2	 33.2	 35		
9	PAT_11_C.fsa	D7S820	 8	 10			
0	PAT_11_C.fsa	CSF1PO	 7	 11			
1	PAT_11_F.fsa	D8S1179	 13	 16			
2	PAT_11_F.fsa	D21S11	 24.2	 33.2			
3	PAT_11_F.fsa	D7S820	 8	 11			
4	PAT_11_F.fsa	CSF1PO	 10	 11			
5	PAT_11_M.fsa	D8S1179	 13	 14			
6	PAT_11_M.fsa	D21S11	 32.2	 33.2			
7	PAT_11_M.fsa	D7S820	 8	 10			
8	PAT_11_M.fsa	CSF1PO	 7	 12			
9	PAT_12_C.fsa	D8S1179	 10	 11			
0	PAT_12_C.fsa	D21S11	 28	 29			
1	PAT_12_C.fsa	D7S820	 10	 12			

Allele Report Settings

Report Style

☐ Allele List

☒ Forensics

☐ Bin Table

☐ Peak Table

☐ Allele Count

☐ Sample Name ☒ File Name

Orientation

☐ Horizontal ☒ Vertical

Options

☒ Extend Diploid Homozygous


☒ Show Allele Name

☐ Show Size (0.1bps)

☐ Show Height

☐ Show Area


☐ Show Score

☒ Show  when no allele call


☐ Show Only Uncertain Alleles

☒ Show Rejected Low Score Alleles

☐ Hide Extra Sample Names



Post Genotyping Applications




The Database search tool allows the user to create a database of profiles, and then search it for possible matches or relatives using Identity by Descent (IBD) to calculate likelihood ratios.

Individual ID : 1
Individual Name: 1
Sample File: PAT_10_C2.fsa

File Name	ID	Name	X/Y	Matched Alleles	Matched Markers	PI/KI
Same-Individual						
PAT_10_C2.fsa	1006		XX	32 32	16 16	3.44E+20
Father/Son						
Mother/Daughter						
Full-Sibs						
PAT_10_C3.fsa	1007		XX	20 32		2.56E+03
PAT_10_C1.fsa	1005		XX	18 32		1.41E+02
Half-Sibs						
PAT_2_F.fsa	682		XY	10 32		1.80E+04
PAT_10_F.fsa	1008		XY	17 32		2.36E+03
PAT_5_F.fsa	691		XY	12 32		7.90E+00
PAT_3_C.fsa	684		XY	13 32		5.77E+00
PAT_1_M.fsa	680		XX	15 32		4.45E+00
PAT_6_F.fsa	694		XY	11 32		3.98E+00
Individ_XYZ	1002		XY	12 32		3.34E+00
PAT_3_F.fsa	685		XY	11 32		3.20E+00
PAT_7_F.fsa	697		XY	12 32		2.06E+00
PAT_8_F.fsa	1022		XY	12 32		1.55E+00

The program will display the most likely relative for a variety of relationships. Likelihood ratios (LR = PI / KI) are displayed for each profile.

The Relationship testing application enables extensive paternity testing using AABB paternity and motherless case calculations. Select the calculation details tab export PI, POP and genotypes for each family.



The built-in Kinship analysis tool allows the user to calculate LR's across three generations.

Kinship Analysis

Single Pair | Multi-Pair

Marker	Individual A			Individual B			Parent/Child (LR)	Full-Sibs (LR)	Half-Sibs (LR)	Uncle-Nephew (LR)	Cousins (LR)	Grandpa
	PAT_10_C3.fsa			PAT_10_C2.fsa								
CSF1PO	11	12		11	13		1.00390	0.75195	1.00195	1.00195	1.00097	1.00195
TPDX	8			8			2.68752	3.39945	1.84376	1.84376	1.42188	1.84376
TH01	7	8		7	8		1.88446	2.72596	1.44223	1.44223	1.22112	1.44223
vWA	16	17		15	18		0.00278	0.25000	0.50000	0.50000	0.75000	0.50000
D16S539	10	12		10	13		2.14998	1.32499	1.57499	1.57499	1.28750	1.57499
D7S820	8	9		9	10		2.30351	1.40176	1.65176	1.65176	1.32588	1.65176
D13S317	10	12		12	14		0.58904	0.54452	0.79452	0.79452	0.89726	0.79452
D5S818	11			11	12		2.14998	1.32499	1.57499	1.57499	1.28750	1.57499
FGA	20	26		22	25		0.00390	0.25000	0.50000	0.50000	0.75000	0.50000
D8S1179	12	13		11	13		1.15181	0.82590	1.07590	1.07590	1.03795	1.07590
D18S51	13	16		13	16		6.34558	18.52269	3.67279	3.67279	2.33640	3.67279
D21S11	27	32.2		27	32.2		7.52494	31.74703	4.26247	4.26247	2.63123	4.26247
D3S1358	15	17		15	17		2.04387	3.28457	1.52194	1.52194	1.26097	1.52194
D2S1338	17	23		17	20		2.51965	1.50983	1.75983	1.75983	1.37991	1.75983
D19S433	14.2	16.2		14	14.2		3.14624	1.82312	2.07312	2.07312	1.53656	2.07312
Product Score:							3.08E-01	2.56E+03	2.03E+02	2.03E+02	3.29E+01	2.03E+02

The preloaded mixture analysis tool can completely deconvolute mixtures with two contributors, and calculate PI, PE, LR, and RMNE

The screenshot displays the Mixture Analysis software interface. On the left, a tree view shows the analysis setup for 'MIX05case2_e' with one contributor. The main table lists markers and their associated data. Below the table, fields for Contributor 1, Contributor 2, Average Major Mx, Cumulative LR, and RMNE are visible. On the right, a 'Trace Data Report' table provides a detailed summary of the analysis results, including PI, PE, LR, and RMNE values for each marker and a cumulative total.

No.	Marker	Major	Minor	Major Mx	Residual	Major Mx
1 -- 1	D8S1179	11,13	11,14	0.87	0.0000	0.95
1 -- 2	D8S1179	11,13	14,14	0.93	0.0085	0.8E
1 -- 3	D8S1179	11,13	13,14	0.88	0.0091	0.7E
2 -- 1	D21S11	28,32,2	30,31	0.83	0.0098	0.74
3 -- 1	D7S820	8,10	9,11	0.87	0.0011	0.9E
4 -- 1	CSF1PO	7,10	12,13	0.90	0.0033	0.8E
5 -- 1	D3S1358	15,15	15,16	0.88	0.0000	
5 -- 2	D3S1358	15,15	16,16	0.94	0.0092	
6 -- 1	TH01	7,9,3	8,10	0.89	0.0092	0.74
7 -- 1	D13S317	12,14	8,9	0.85	0.0032	0.84
8 -- 1	D16S539	10,11	9,12	0.85	0.0010	0.9E
9 -- 1	D2S1338	17,21	16,24	0.90	0.0023	0.94
10 -- 1	D19S433	13,13	13,14	0.84	0.0003	
10 -- 2	D19S433	13,13	12,14	0.82	0.0053	
10 -- 3	D19S433	13,13	14,14	0.92	0.0054	
11 -- 1	vWA	15,15	16,19	0.89	0.0006	
12 -- 1	TPDX	9,10	8,11	0.86	0.0062	0.7E
13 -- 1	D18S51	17,18	12,15	0.84	0.0012	1.0C

No.	Marker	Mixture	PI	PE	Contributor 1	Contributor 2	LR
1	D8S1179	11,13,14	0.31586	0.68414	11,14	11,13	51.685
2	D21S11	28,30,31,32,2	0.32685	0.67315	30,31	28,32,2	33.365
3	D7S820	8,9,10,11	0.77414	0.22586	9,11	8,10	6.381
4	CSF1PO	7,10,12,13	0.41470	0.58530	12,13	7,10	37.064
5	D3S1358	15,16	0.40653	0.59347	15,16	15	10.941
6	TH01	7,8,9,3,10	0.53097	0.46903	8,10	7,9,3	11.361
7	D13S317	8,9,12,14	0.27584	0.72416	8,9	12,14	33.775
8	D16S539	9,10,11,12	0.68160	0.31840	9,12	10,11	13.529
9	D2S1338	16,17,21,24	0.14541	0.85459	16,24	17,21	35.002
10	D19S433	12,13,14	0.34027	0.65973	12,14	13	16.508
11	vWA	15,16,19	0.24614	0.75386	16,19	15	28.890
12	TPDX	8,9,10,11	0.73706	0.26294	8,11	9,10	31.457
13	D18S51	12,15,17,18	0.26381	0.73619	12,15	17,18	26.882
14	AMEL	X,Y	--	--	X	X,Y	--
15	D5S818	8,11,12,13	0.76054	0.23946	11,12	8,13	43.294
16	FGA	20,23,24	0.12168	0.87832	23,24	20,24	72.871
		Cumulative:	2.95E-07	1-2.95E-07			1.03E+21

Contributor 1: MIX05case2_victim (Minor) ☐ Contested

Contributor 2: Individ_XYZ (Major) ☒ Contested

Average Major Mx: 0.870

Cumulative LR: 1.03E+21 RMNE: 1-2.95E-07

Comment:

If no contributors are found in the current project, the user may search GeneMarkerHID's built in database.

More Information

Thank you for your interest in GeneMarker HID

STR Human Identity Software

Have more questions? Please email tech_support@softgenetics.com or call 814-237-9340 to...

- Obtain more information about GeneMarkerHID
- Request a Quotation
- Schedule a **Free**, online training session with a SoftGenetics expert

Finally, please note that the **User Manual** is always available by navigating to **Help → Help**.

